

Any equipment we provide is on loan and as such should be cared for in accordance with the instructions and re- turned to the Integrated Community Equipment Loans Store (ICELS) by telephoning the number below.

If you discover any damage, defects or problems operating the equipment then you must report this immediately to

British Red Cross –

Telephone:
0345 127 2911

Email:
customercare@redcross.org.uk



ICELS

INTEGRATED COMMUNITY EQUIPMENT LOAN SERVICE



What to do if you are concerned

Please talk to any health and social care professional involved in your care or contact the number below

Nottinghamshire County Council
Customer Service Centre [0300 500 80 80](tel:03005008080)

Safe Use of Bed Levers in the Community

ICELS

INTEGRATED COMMUNITY EQUIPMENT LOAN SERVICE



Tel: 0345 127 2911

Safe Use of Bed Levers

This bed lever has been provided to help you when turning in bed, and when getting out of the bed.

A risk assessment has been carried out by your Health and Social Care professional to assess if this equipment will be of benefit to you.

You will be given a copy of your bed lever risk assessment. It is important that you

keep this information safe and available so that every person involved in your care understands your needs.


This bed lever has been assessed as compatible with your current bed. It has been fitted to comply with both the manufacturer's instructions and national guidance on the safe use of bed levers. It is important that you do not change the position of the bed lever without consulting a health and social care professional with specific training in bed levers.

Most manufacturer's state that a minimum gap of 318mm (12.5 inches) is essential between the headboard and the start of the rail on the bed lever.

You are advised that it is your responsibility / carers responsibility to regularly check the piece of equipment.

Only use the bed lever for its intended purpose.

To maximise your safety, do not move it from the original position.



Regular Checks –what to look for

All parts of the equipment are in safe working order

Check for defects. If you notice anything has become loose, seems to be cracked, worn, or missing, then bed rail may not be safe to use.

Cleaning

All parts of the equipment can be cleaned with a mild detergent using a cloth. Avoid letting water run inside the tubes.

Using Bed Levers

It is important that you and your carer/s monitor the use of the equipment and report any of the following changes-

- Changes in the ability of the person using the equipment which may impact on the equipment being used safely.
- Any change of bed. Your bed lever has been assessed as safe with your current bed. If the bed is changed, the bed lever may no longer be safe. You may need an alternative design of bed lever.
- If the equipment moves out of position, resulting in a different gap between the bed lever /headboard/ mattress
- If there are concerns about any other gaps that could present as an entrapment risk.
- If the bed lever needs refitting, for example, if you need to move the bed into a different room.
- If you are getting caught against the bed lever and can't independently move away from it.
- If you have trapped any part of your body in, or at the side of the bed lever.
- If you develop new involuntary movements in your arms and legs.
- If you are confused or more confused than usual